



Air conditioning

Description

The DLB 2000, well established in practice, ensures fully automated regulation of rel. humidity with top stability. Input of the required rel. humidity and setting of the temperature takes place using a touchscreen display. A regulated flow of dry air has to be supplied at the DLB 2000 inlet. The precision of the adjustable humidity is guaranteed with a mass flow controller which precisely adjusts the air flow. In the DLB 2000 the air is precisely regulated to the specified value of rel. humidity. The DLB 2000 has a large water reservoir, which can be heated up to 50 °C, sufficient for up to 8 hours of operation. This function guarantees extremely stable adjustment of rel. humidity and temperature. Additional fields of application for the DLB 2000 are laboratory tasks with high requirements regarding the stability of rel. humidity and air temperature such as adsorption experiments on filter media.

Benefits

- Adjustable relative humidity from 20 to 80 %
- Highest stability of relative humidity 0.2 %
- Easy handling by touch screen
- External control analogue 0 - 5 V
- Huge water supply tank
- Easy refilling of water
- Reliable function
- Low maintenance
- Reduces your operating costs!

Datasheet

<i>Parameter</i>	<i>Description</i>
Interfaces	Analog (0 – 5 V (set point), 0 – 5 V (measured humidity))
Volume flow	10 – 60 NI/min (others on request)
Power supply	115/230 V, 50/60 Hz
Dimensions	280 • 320 • 460 mm (H • W • D)
Weight	approx. 16 kg
Maximum counter pressure	0.5 bar _g overpressure
Luftfeuchteregelung	20 – 80 % +/- 0.2 % at 23 °C
Water supply	1 l deionised water

Applications

- Accurate regulation of rel. humidity for test rigs, e.g.
 - filter testing
 - active carbon filter testing
 - bioaerosol measuring device testinn
 - dust measuring device testing
 - etc.
- Conditioning of the dispersion air for the RBG 1000 and RBG 2000 dust dispersers (the instruction manual of the dust dispersers has to be followed).
- Laboratory experiments with special requirements for the stability of rel. humidity.

Palas GmbH
Partikel- und Lasermesstechnik
Greschbachstrasse 3 b
76229 Karlsruhe
Germany

Managing Partner:
Dr.-Ing. Maximilian Weiß
Commercial Register:
register court: Mannheim
company registration number: HRB 103813
USt-Id: DE143585902



Contact: E-Mail: mail@palas.de Internet: www.palas.de Tel: +49 (0)721 96213-0 Fax: +49 (0)721 96213-33