






Microparticles with a mean diameter of 40 µm

Description

MP 40 are microparticles with a mean diameter of 40,06 µm. The distribution of this dust is monodisperse, the geometric standard deviation of the distribution is only 3,56 µm. Due to this small distribution these dusts are used for scientific analysis and can also be used for the calibration of measuring devices. The particles consist of PMMA (Polymethylmethacrylat).





Benefits

- Feststoffpartikel mit sehr enger monodisperser Verteilung
- Hohe Standzeit des erzeugten Aerosols

Datasheet

<i>Parameter</i>	<i>Description</i>
Particle material	PMMA
Particle diameter (maximum)	75 µm
Particle diameter (minimum)	0.5 µm
Particle diameter (mean)	40.06 µm
Geometric standard deviation (volume)	3.56 µm

Applications

Palas GmbH
Partikel- und Lasermesstechnik
Greschbachstrasse 3 b
76229 Karlsruhe
Germany

Managing Partner:
Dr.-Ing. Maximilian Weiß
Commercial Register:
register court: Mannheim
company registration number: HRB 103813
USt-Id: DE143585902



Contact: E-Mail: mail@palas.de Internet: www.palas.de Tel: +49 (0)721 96213-0 Fax: +49 (0)721 96213-33